

# CERAGYR

THE MOBILE-BEARING KNEE PROSTHESIS



THE ULTIMATE CONGRUENCY

 **CERAVER**  
EXPERIENCE - INNOVATION

## CERAGYR: THE MOBILE-BEARING KNEE PROSTHESIS

## BENEFITS

Excellent mobility, particularly in rotation and flexion

Full stability in flexion/extension

Optimal resistance to creep and wear

Stable, unconstrained patello-femoral articulation

Full interchangeability of component sizes




## THE ULTIMATE CONGRUENCY


## DESIGN

**Femoral component with**


- Anatomic trochlear depth and alignment, facilitating patellar contact and tracking.
- Distal and posterior condyles which have a spherical curvature radius.

**"Spherical" patellar button**

allows improved contact with trochlea and ensures excellent patello-femoral stability.

**Freely rotating mobile bearing**

and anterior-posterior motion reduce stresses at the various interfaces.

**Tibial tray anchorage to base-plate**

eliminates all risk of dislocation.



CERAGYR: THE MOBILE-BEARING KNEE PROSTHESIS

1997 - 1999

Research studies


SPECIAL FEATURES

ULTIMATE CONGRUENCY

The CERAGYR prosthesis features "spherical" femoral-tibial contact surfaces associated with a mobile bearing. This allows perfect congruency of bearing surfaces, accompanied by low stress levels, the ideal combination for long-term life expectancy of the polyethylene bearing and the bond.

◆ "Spherical" femoral-tibial contact surfaces

An optimal spherical 27 mm radius ensures stable anterior-posterior and medio-lateral articulation, with perfect conformity of trochlear groove and condyles.




◆ Maximisation of contact surface areas

This results in significantly reduced contact pressure (3 Mpa, which is substantially below the polyethylene elastic limit) and reduced creep and wear

Prosthesis	A	B	C	CERAGYR
Flexion	Surface (mm <sup>2</sup> )	Surface (mm <sup>2</sup> )	Surface (mm <sup>2</sup> )	Surface (mm <sup>2</sup> )
0°	731	1142	457	1161
20°	843	786	264	1140
40°	829	742	182	1136
60°	748	784	177	1145
80°	486	735	153	800

◆ Complete interchangeability of component sizes

The use of a single radius and a fixed interaxial distance between spheres allows complete interchangeability of femoral components and tibial trays. Ensuring a better match to patient anatomy while reducing surgery times.




MINIMISED STRESSES

◆ Free rotation of the tibia

The freely rotating mobile bearing eliminates torsional stresses, encouraging bone ingrowth without the need for precise alignment of the tibial and femoral components.

◆ Antero-posterior glide of the tibia preserved

The bearing allows antero-posterior glide of 5 mm so that its position relative to the tibial base-plate in the antero-posterior plane is determined primarily by the patellar reaction force and the position of the joint line, eliminating purely rotational stresses.



■ AN ORIGINAL DESIGN CONCEPT, MASTERFULLY IMPLEMENTED

1999 - 2004

Clinical validation of the implant and associated surgical technique

## MOBILITY WITHOUT RISK OF DISLOCATION

### ◆ Mobility in flexion


The condyles' rotation centre is located behind the bearing center. The resultant reduction of the posterior lip height allows a greater degree of flexion.

### ◆ Medial tibial spine

At high flexion any risk of tray dislocation under the femoral component is limited.

### ◆ Tibial tray anchored to base-plate

The tibial tray is equipped with a ratchet mechanism on the tibial peg, which eliminates the risk of tray dislocation.



## LONG TERM FIXATION

### ◆ Rectangular cross-section keel with fins

The keel offers an excellent anchorage for the implant. The fins, with a posterior alignment of 30°, ensure that the base-plate remains securely in place under rotation, avoiding the risk of tibia fracture.

### ◆ Modular tibial base-plate

The CERAGYR base-plate is compatible with a variety of tibial extension options, both straight-stemmed and eccentric, for use in primary and revision procedures.



The femoral component and the tibial base-plate of the CERAGYR prosthesis are both available in cemented or uncemented versions. Their shape and dimensions are identical, allowing the choice to be deferred until the end of the procedure.

### ◆ Cemented version

The implant undersurfaces incorporate voids capable of accommodating cement under pressure during compaction, so ensuring that the implant is securely anchored.

### ◆ Uncemented version

The implant undersurfaces are coated with porous beads whose pore size (250µm) is ideal for the bone ingrowth that ensures optimal secondary implant fixation.



## CERAGYR: THE MOBILE-BEARING KNEE PROSTHESIS

## SPECIAL FEATURES

## THE PATELLO-FEMORAL PROSTHESIS, THE REFERENCE

A design tried and tested since 1979 with the **HERMES** range.

The design of the patello-femoral articulation in the **CERAGYR** prosthesis ensures optimal contact and centring of the patella within the trochlea over the full range of motion and minimal stresses at the various interfaces.

◆ **Near-anatomic trochlea**

The trochlear groove is deep with a  $6^\circ$  alignment upward and outward. This design results in reduced stress within the natural or prosthetic patella.

◆ **Single radius of curvature for patellar button and trochlea**



This geometry ensures optimal contact whatever the component size, allowing complete interchangeability.

◆ **Elevated outer face**

This helps to reduce the risk of patellar displacement and limits peak stresses within the patellar button.

◆ **Optimisation of patellar contact with trochlea**

At full extension the trochlea is raised sufficiently to ensure contact with the patella. The trochlear groove is inclined as far as the top of the anatomic intercondylar notch to ensure optimal contact even at maximum flexion



## 2-YEAR RESEARCH STUDIES AND 6-YEAR CLINICAL VALIDATION

◆ **Effect of wear**

Simulation studies have found polyethylene wears less than or equal to 0.05 mm/year, with contact surfaces greater than 1,000 mm<sup>2</sup>. The polyethylene ratchet mechanism shows no signs of deterioration after more than 10 million cycles. These results guarantee the reliability of the solutions developed. Initial clinical results confirm in-vitro findings.

◆ **Clinical evaluation 1999–2004**

Radiological studies have shown that in over 50% of cases the tibial tray does not impinge on the spine until flexion has reached 90°. This observation shows the importance of ensuring antero-posterior motion is preserved even by a maximum of 5 mm, to reduce the shear stresses at the bonding interfaces.





### PRODUCT REFERENCES

#### PATELLAR BUTTON UHMWPE\* (ISO 5834/1-2)

Size	Ø 31	Ø 34	Ø 36	Ø 39
Thickness	8 mm	8 mm	9 mm	11 mm
Ref	<b>5090</b>	<b>5091</b>	<b>5092</b>	<b>5093</b>



#### FEMORAL COMPONENT - cobalt chromium alloy (ISO 5832/4)


Size	1	2	3	4	5	
DIM.	AP (mm)/Box	53/35	57/39	61/43	65/47	69/51
	ML (mm)	60	64	68	72	77

#### CEMENTED

Ref	Right	6201	6202	6203	6204	6205
	Left	6101	6102	6103	6104	6105

#### UNCEMENTED

Ref	Right	6611	6612	6613	6614	6615
	Left	6111	6112	6113	6114	6115



#### TIBIAL TRAY - UHMWPE\* (ISO 5834/1-2)

Size	1	2	3	4 - 5	
Ref	Th. 6 mm	6311	6312	6313	6314
	Th. 8 mm	6321	6322	6323	6324
	Th. 10 mm	6331	6332	6333	6334
	Th. 12 mm	6341	6342	6343	6344
	Th. 14 mm	6351	6352	6353	6354



#### TIBIAL BASE-PLATE - cobalt chromium alloy (ISO 5832/4)


Size	1	2	3	4	5	
DIM.	AP (mm)	41	43	47	51	55
	ML (mm)	64	68	73	78	84

#### CEMENTED

Ref.	6401	6402	6403	6404	6405
------	------	------	------	------	------


#### UNCEMENTED

Ref.	6501	6502	6503	6504	6505
------	------	------	------	------	------



#### TIBIAL EXTENSION (length 50 mm) - titanium alloy (ISO 5832/3)

Size	Ø10	Ø11	Ø12	Ø13	Ø14	Ø15	Ø16	Ø17	Ø18	Ø19
Ref.	<b>5080</b>	<b>5081</b>	<b>5082</b>	<b>5083</b>	<b>5084</b>	<b>5085</b>	<b>5086</b>	<b>5087</b>	<b>5088</b>	<b>5089</b>



\* All UHMWPE implants are manufactured from compression moulded sheeting and sterilised by gamma irradiation in an argon gas environment (minimum dose 25 kGy).

## THE CERAVER KNEE RANGE

### HERMES, a total range

HERMES FP



The patella-femoral

HERMES UNI



The unicompartmental

HERMES 1C-2C



Both cruciates or posterior cruciate retaining

HERMES PS



The posterior-stabilized

HERMES REVISION



The revision posterior-stabilized

### CERAGYR

Ultimate congruency



The mobile bearing knee

### CERAVISON

Computer-assisted surgery system



PARIS NORD 2 - 69, RUE DE LA BELLE ÉTOILE - B.P. 54263 GONESSE  
95957 ROISSY CDG CEDEX - FRANCE

TEL. : (+33) 1 48 63 88 63 - FAX : (+33) 1 48 63 88 99

WWW.CERAVER.COM - CONTACT@CERAVER.FR

CERAVER OSTEAL AMERICA INC.  
2277 GUENETTE  
SAINT LAURENT QUEBEC  
H4R 2E9 CANADA  
TEL : 19 1 514 856 24 35  
FAX : 19 1 514 856 24 44

CERAVER DEUTSCHLAND  
THEODOR HEUSS STRASSE 52  
61118 BAD VILBEL  
ALLEMAGNE  
TEL : 00 49 6101 829844  
FAX : 00 49 6101 829846

CERAVER IBERICA  
BELTRAN BAGUENA 5-10<sup>a</sup>-10A  
NUEVO CENTRO - 46009 VALENCIA  
ESPAGNE  
TEL : 00 34 9 6 348 10 19  
FAX : 00 34 9 6 348 10 02

CERAVER ITALIA  
VIA DELL'INDUSTRIA, 2 - BARGELLINO  
40012 CALDERARA DI RENO  
BOLOGNA - ITALIE  
TEL : 00 39 051 72 88 76  
FAX : 00 39 051 72 78 03